



## **MST Simulated Vomit: Usage & Instructions**

Made from the same silicone rubber used in movie productions, this simulated vomit is safe to use in any training scenario. Durable and reusable, it will last through thousands of trainings.

### ***Usage:***

MST Simulated Vomit is a standalone product and can be used in any training scenario. Lay MST Simulated Vomit on any flat surface to create the illusion of fresh vomit. MST Simulated Vomit is compatible with MST Vomit Odor.

### ***Cleaning:***

MST Simulated Vomit may be rinsed in soapy warm water to remove residue, simulated blood, and other contaminants.

### ***Warnings:***

Do not stretch MST Simulated Vomit with excessive force as rips may occur. Avoid stepping on MST Simulated Vomit as degradation may occur.

### ***Related Products:***

- MST Vomit Odor